AB 2886 Guidance Letter No. 016: 8260FAB: Fuel Oxygenates and Additives plus BTEX Minimal List

State Water Resources Control Board (SWRCB)

Letter No. 016

Subject: 8260FAB; GC/MS Analysis of Fuel Oxygenates and Additives

plus BTEX

Date: April 26, 2002

Overview:

Within the LUST Cleanup program, many soil and water samples are analyzed for VOCs and oxygenates by GC/MS USEPA Method 8260B. Electronic EDF data for these analyses are reported to the SWRCB GeoTracker system as required by AB2886.

For data comparability purposes the SWRCB recommends the following minimal analyte list, spiking compounds, surrogate compounds and reporting limits for this method. This analyte list is includes of a minimal list of fuel oxygenates and BTEX compounds. USEPA Method 8260B is recommended for analysis of Tert-butyl Alcohol.

Special Conditions:

This applies to all sample matrices.

Areas of Impact:

Field(s): ANMCODE, PARLABEL

Entry: ANMCODE = "8260FAB" (Volatile Organic Compounds by GC/MS Fuel

Additives plus BTEX)

PARLABEL = various; refer to list below

Policv:

a) Quality Control Requirements

Preparation/analytical batch should include:

- one laboratory method blank
- one matrix spike
- one matrix spike duplicate or one matrix duplicate if target compounds are present
- one blank spike

Control limits for the blank spike/blank spike duplicates are 70% - 130% per SW-846.

AB 2886 Guidance Letter No. 016: 8260FAB: Fuel Oxygenates and Additives plus BTEX Minimal List

b) Analytes

"8260FAB" Spiking Compounds¹

| PARLABEL | Description | | | | |
|------------------------------------|---|--|--|--|--|
| Halogenated Compounds ¹ | | | | | |
| DCA12 | 1,2-Dichloroethane ¹ | | | | |
| ETHANOL | Ethanol (EtOH) ¹ | | | | |
| TBA | tert-Butyl alcohol (TBA) ¹ | | | | |
| Aromatic Compounds ¹ | | | | | |
| BZ | Benzene ^{1,2} | | | | |
| BZME | Toluene ^{1,2} | | | | |
| Oxygenate Compounds ¹ | | | | | |
| MTBE | Methyl-tert-butyl ether (MTBE) ¹ | | | | |

¹ From Region 4

"8260FAB" Analytes and Detection Limits

| | | | | Water (ug/L) | | | Soil (ug/kg) | |
|----------|------------------------------------|-------|-------|--------------|------|-----|--------------|------|
| PARLABEL | Description | Order | Surr. | MDL | MRL | MCL | MDL | MRL |
| MTBE | Methyl-tert-butyl ether (MTBE) | 1 | | | 2 | | | 5 |
| ETBE | Ethyl tert-butyl ether (ETBE) | 2 | | | 2 | | | 5 |
| TAME | tert-Amyl methyl ether (TAME) | 3 | | | 2 | | | 5 |
| DIPE | Di-isopropyl ether (DIPE) | 4 | | | 2 | | | 5 |
| TBA | tert-Butyl alcohol (TBA) | 5 | | | 10 | | | 20 |
| DCA12 | 1,2-Dichloroethane | 6 | | | 0.5 | 0.5 | | 1.0 |
| EDB | 1,2-Dibromo ethane | 7 | | | 1 | | | 2 |
| ETHANOL | Ethanol (EtOH) | 8 | | | 1000 | | | 1000 |
| BZ | Benzene | 9 | | 0.04 | 0.5 | 1.0 | | 1.0 |
| BZME | Toluene | 10 | | 0.11 | 1 | | | 2 |
| EBZ | Ethylbenzene | 11 | | 0.06 | 1 | | | 2 |
| XYLENES | Xylenes | 12 | | | 3 | | | 6 |
| DBFM | Dibromofluoromethane ¹ | | 4 | | | | | |
| DCA12D4 | 1,2-Dichloroethane-d4 ¹ | | 4 | | | | | |
| BR4FBZ | 4-Bromofluorobenzene ¹ | | 4 | | | | | |
| BZMED8 | Toluene-d8 ¹ | | 4 | | | | | |

 $MDL-Method\ Detection\ Limit\ (from\ SW-846\ for\ wide-bore\ capillary\ columns)$

² Minimum spiking requirements recommended in SW-846

MRL – Method Reporting Limit (from Region 4)

 $MCL-Maximum\ Contamination\ Limit\ (from\ DHS\ California\ Drinking\ Water\ Quality\ Database\ 01/08/2002)$

¹ SW-846 recommended surrogate