

SECTION 9
RADIOACTIVE STRONTIUM IN DRINKING WATER
Method 905.0

1. Scope and Application

- 1.1 This method covers the measurement of total strontium and soluble strontium-89 and strontium-90 in drinking water. Some naturally insoluble (and probably suspended) forms of strontium-89 and strontium-90 would also be measured by this method when samples of such drinking water supplies are acid-preserved before analysis.
- 1.2 The Drinking Water Regulations under the Safe Drinking Water Act set maximum contaminant concentrations for radionuclides in drinking water based on a 2 liter per day drinking water intake using the 168 hour data listed in Handbook 69, National Bureau of Standards. The maximum contaminant concentration for strontium-89 and strontium-90 are 80 pCi/l and 8 pCi/l, respectively, the critical organ being bone marrow. If other radionuclides are also present in the drinking water, the sum of their annual dose equivalent must not exceed 4 mrem per year. The Regulations also give a required sensitivity of measurement which is defined in terms of a detection limit. The required detection limits given for strontium-89 and strontium-90 are 10 pCi/l and 2 pCi/l, respectively. Appendix C has equations for calculating the counting time necessary to meet the required detection limit.

2. Summary of Method

- 2.1 Stable strontium carrier is added to the drinking water sample and strontium-89 and strontium-90 are precipitated from the solution as insoluble carbonates. Interferences from calcium and some radionuclides are removed by one or more precipitations of the strontium carrier as strontium nitrate. Barium and radium are removed as the chromate. The yttrium-90 daughter of strontium-90 is removed by a hydroxide precipitation step and the separated combined strontium-89 and strontium-90 are counted for beta particle activity. The counting result; immediately ascertained, represents the total strontium activity (strontium-90 + strontium-89) plus an insignificant fraction of the yttrium-90 that has grown into the separated strontium-90. The yttrium-90 daughter grows in again and is then separated with stable yttrium carrier as hydroxide and finally precipitated as oxalate and beta counted. The strontium-90 concentration is determined by the yttrium-90 activity and the strontium-89 concentration is then determined by difference.
- 2.2 Counting efficiency data must be obtained with standard strontium-89, strontium-90, and yttrium-90 activities. These data are used to make corrections since strontium-89, strontium-90, and yttrium-90 emit beta particles with different energies.

3. Sample Handling and Preservation

3.1 It is recommended that samples be preserved with acid at the time of collection. For preservation, sufficient acid should be added to make the sample $\text{pH} < 2$.

3.2 The Drinking Water Regulations allow for the option of quarterly compositing for an annual analysis or averaging the analyses of four quarterly samples. It is especially recommended to preserve composited samples.

3.3 It is recommended that no less than one liter size samples be collected for analysis.

4. Interferences

4.1 Radioactive barium and radium will be carried down with radioactive strontium as carbonate. This method includes steps to separate strontium from barium and radium.

4.2 Samples that naturally contain significant amounts of stable strontium will cause errors in the recovery of the added strontium carrier. Blank samples to which no strontium carrier is added should be run to determine natural strontium content. Hard waters contain calcium which precipitates with the strontium in the initial carbonate precipitation. If not separated, the calcium will cause errors in the recovery of the strontium carrier. Repeated precipitations with 16N HNO_3 (conc.) will eliminate this interference.

5. Apparatus - See Appendix D for details and specifications.

5.1 Low background beta counting system (< 3 cpm background on the beta voltage plateau is recommended).

5.2 Centrifuge and 50 ml centrifuge tubes

5.3 Drying oven

5.4 Hot water bath

5.5 Electric hot plate

5.6 Analytical balance

5.7 pH meter

5.8 Desiccator, aluminum and/or glass

5.9 Stainless steel planchets, 2-inch diameter by 1/4-inch deep

- 5. 10 Sintered-glass (fine). crucibles
- 5.11 Plastic ring and disc mounts
- 5.12 Mylar film
- 5.13 Teflon filter holder
- 5. 14 Drying lamps
- 5.15 Glassware

6. Reagents

Distilled or deionized water is to be used, and all chemicals should be of "reagent-grade" or equivalent whenever they are commercially available.

- 6.2 Strontium carrier (10 mg/ml): Dissolve 24.16g $\text{Sr}(\text{NO}_3)_2$ in water and dilute to 1 liter in a volumetric flask with water. Mix thoroughly.

Standardization: (In triplicate).

Carefully pipet 10.0-ml portions of the strontium carrier solution into separate 50-ml centrifuge tubes. Add 1 ml 6N NaOH and heat in a water bath. Slowly, and with stirring, add 15 ml of 2N Na_2CO_3 solution (see sodium carbonate solution below) and continue digesting for 15 to 20 minutes. Allow to cool and filter the SrCO_3 precipitate through a tared sintered-glass (fine) crucible. Wash the precipitate and the crucible walls with three 5-ml portions of distilled water adjusted to pH 8 with 6N NH_4OH , and with three 5-ml portions of acetone. Dry the crucible for 30 minutes in a 105°C oven. Cool the crucibles in a desiccator and weigh.

$$\text{Strontium, mg/ml} = \frac{\text{mg of SrCO}_3 * 0.5935}{10\text{ml}}$$

- 6.3 Yttrium carrier (10 mg/ml): Dissolve 43g $\text{Y}(\text{NO}_3)_3 \cdot 6\text{H}_2\text{O}$ in water plus 5 ml 16N HNO_3 (conc.) and dilute to 1 liter in a volumetric flask with water. Mix thoroughly

Standardization: (In triplicate).

Carefully pipet 10.0-ml portions of the yttrium carrier solution into separate 50-ml centrifuge tubes. Add 30 ml saturated $(\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$ to each centrifuge tube and stir. Digest in a hot water bath (near boiling) for

30 minutes. Cool in an ice bath. Filter the precipitate onto a Whatman #42 filter paper, then ignite in a tared crucible at 800°C for 1 hour to convert the oxalate to the oxide. Cool and weigh the crucible and calculate the yttrium concentration from the following equations.

$$\frac{Y}{Y_2O_3} = \frac{2 * 88.92}{225.84} = 0.7875$$

$$\text{yttrium, mg/ml} = \frac{\text{mg of } Y_2O_3 * 0.7875}{10\text{ml}}$$

- 6.4 Acetic acid, 5.8N: Mix 1 volume 17.4N CH₃COOH (glacial) with 2 volumes of water.
- 6.5 Acetone, (CH₃)₂CO: anhydrous.
- 6.6 Ammonium acetate buffer: Dissolve 154g NH₄C₂H₃O₂ in 800 ml of water, add 57 ml 17.4N CH₃COOH (glacial), adjust the mixture to pH 5.5 using CH₃COOH or NH₄OH. Dilute to 1 liter.
- 6.7 Ammonium hydroxide, 15N: NH₄OH (conc.), sp. gr. 0.90, 56.6%.
- 6.8 Ammonium hydroxide, 6N: Mix 2 volumes 15N NH₄OH (conc.) with 3 volumes of water.
- 6.9 Ammonium hydroxide, 0.1N: Mix 1 volume 15N NH₄OH (conc.) with 150 volumes of water.
- 6.10 Ammonium oxalate, saturated: Into 100 ml boiling water, dissolve 10g (NH₄)₂C₂O₄·H₂O. Cool.
- 6.11 Barium carrier, (10 mg/ml): Dissolve 19.0g Ba(NO₃)₂ in water and dilute to 1 liter with water.
- 6.12 Hydrochloric acid, 6N: Mix 1 volume 12N HCl (conc.) with 1 volume of water.
- 6.13 Methyl red indicator, 0.1%: Dissolve 0.1g of methyl red in 100 ml ethanol.
- 6.14 Nitric acid, 16N: HNO₃ (conc.) sp. gr. 1.42, 70.4%.
- 6.15 Nitric acid, 6N: Mix 3 volumes 16N HNO₃ (conc.) with 5 volumes of water.
- 6.16 Nitric acid, 1N: Mix 1 volume 6N HNO₃ with 5 volumes of water.
- 6.17 Phenolphthalein indicator, 1%: Dissolve 1g phenolphthalein in 50 ml ethanol and

add 50 ml water.

- 6.18 Sodium carbonate, 2N: Dissolve 124g $\text{Na}_2\text{CO}_3 \cdot \text{H}_2\text{O}$ (or 106g Na_2CO_3) in water and dilute to 1 liter with water.
- 6.19 Sodium chromate, 0.5M: Dissolve 117g $\text{Na}_2\text{CrO}_4 \cdot 4\text{H}_2\text{O}$ in water and dilute to 1 liter with water.
- 6.20 Sodium hydroxide, 6N: Dissolve 240g NaOH in water and dilute to 1 liter with water.
- 6.21 Wetting agent solution: e.g. Photo-Flo, Eastman Kodak Co.

7. Calibrations

- 7.1 Counting Efficiencies - Separate counting efficiencies should be determined for strontium-89 and strontium-90 using known amounts of the respective radioactive standards and 20.0 mg of strontium carrier, precipitated as carbonate and counted. A strontium-90 precipitate is prepared after separation of the yttrium-90 daughter by the following procedure. Add a known amount of strontium-90 standard, in the order of 1000 disintegrations per minute (dpm), and 20 mg of strontium carrier to a 50-ml centrifuge tube, add 20 ml of water and proceed as in steps 8.9 through 8.11. Then for the yttrium-90 counting efficiency, continue with steps 8.12 through 8.16.
- 7.2 Sources of supply: For strontium-90-yttrium-90, the National Bureau of Standards, Washington, DC offers a standard solution (SRM 4234) as listed in their latest catalog #260.

For strontium-89. Amersham Radiochemicals, Arlington Heights, Illinois, offers a standardized aqueous solution essentially free from strontium-90. This item is listed as SMZ.64 in their latest catalog.

Standard sources of strontium-89 and strontium-90 are also available from the Quality Assurance Division, U.S. Environmental Protection, EMSL-Las Vegas.

8. Procedure

- 8.1 Transfer 1-liter water sample aliquots to 2-liter beakers. Add 2.0 ml each of strontium and barium carrier solutions to each sample and blank beakers. Heat the samples to boiling and add 6N NaOH while stirring, to the phenolphthalein end point (red color), and add 50 ml 2N Na_2CO_3 solution. Continue heating to near boiling for 1 hour with occasional stirring. Then set the beakers aside for at least 2 hours, allowing the carbonate precipitate to settle.
- 8.2 Decant most of the clear supernate and discard it. With the remainder of the

supernate and necessary water washes (adjusted to pH 8 with 6N NH₄OH), quantitatively transfer the precipitate to a 50-ml centrifuge tube. Centrifuge and discard the supernate. This precipitate will contain the strontium and barium carriers.

- 8.3 Dissolve the precipitate by the dropwise addition of 4 ml 16N HNO₃.
- 8.4 Add 20 ml 16N HNO₃ to the centrifuge tube, cool in an ice bath and stir. Centrifuge and discard the supernate which will contain a significant fraction of the calcium present in the sample.
- 8.5 Add 20 ml 16N HNO₃ to the centrifuge tube, cool in an ice bath and stir. Centrifuge and discard supernate.

Note: If drinking water samples contain much calcium (hardness), it will be necessary to repeat step 8.5.

- 8.6 Dissolve the strontium and barium nitrate precipitate in 25 ml water, add 2 drops methyl red indicator, neutralize to yellow color with 6N NH₄OH, then adjust the pH back to red color by adding 5.8N CH₃COOH dropwise.
- 8.7 Add 5 ml ammonium acetate buffer solution, and heat in a hot water bath. Add, with stirring, 2 ml 0.5M Na₂CrO₄ and digest in the hot water bath for 15 minutes. Cool the reaction mixture and centrifuge. Transfer the supernate to a clean 50-ml centrifuge tube, and discard the barium chromate residue.

Note: This residue can be saved if radioactive barium, radium, or lead analysis is desired.

- 8.8 To the buffered chromate supernate add 2 ml 15N NH₄OH and heat in a hot water bath. Add 5 ml 2N Na₂CO₃ solution and digest for 15 minutes. Cool, centrifuge, and discard the supernate.

Note: In the next step, the strontium-89 and strontium-90 are separated from yttrium-90 with a yttrium carrier scavenge to start a specific ingrowth period and to get a separate radiostrontium count in the following steps.

- 8.9 Add a few drops 16N HNO₃ to the carbonate precipitate, then add 25 ml water and 1 ml yttrium carrier. Add 1 drop of wetting agent solution (such as "Photo-Flo," an Eastman Kodak Company film processing product) and 5 ml 15N NH₄OH. Heat in a hot water bath for 15 minutes with occasional stirring. Centrifuge and transfer the supernate to a clean 50-ml centrifuge tube. Wash the yttrium hydroxide precipitate with 5 ml water, centrifuge and add this wash to the supernate. Note the time of this yttrium hydroxide precipitation which marks the beginning of the yttrium-90 ingrowth period. From this point on it is important to proceed without delay to the final separation and count of the strontium-89 and

strontium-90 activity to minimize ingrowth of yttrium-90.

Note: Concentrated NH_4OH sometimes contains CO_2 in solution which will cause precipitation of some of the strontium carrier in this step. If low carrier recoveries are obtained in step 8.11, then for subsequent strontium analyses, anhydrous NH_3 gas may be substituted for concentrated NH_4OH in step 8.9 by bubbling NH_3 gas in the sample solution until the phenolphthalein end point is reached, and then 5 minutes more. The same precaution might be taken in step 8.14 to prevent carrydown of the strontium-90 as the carbonate precipitate in that step.

- 8.10 Add 5 ml 2N Na_2CO_3 to the supernate from step 8.9, heat in a hot water bath for about 10 minutes, centrifuge and discard the supernate.
- 8.11 Slurry the strontium carbonate precipitate with a few ml water and transfer quantitatively to a tared glass fiber filter. Wash the precipitate with three 10-ml portions of water adjusted to pH 8 with NH_4OH , then with three 10-ml portions of acetone. During filtration and washes of the strontium carbonate, minimize the time of air flow through the filter to avoid collection of radon daughters. Dry the filter at 105°C for 10 minutes, then weigh, mount and count (within 2 hours). This count gives the total of strontium-89 and strontium-90 activities, plus the ingrown yttrium-90. Note the time of this count as it must be corrected for yttrium-90 ingrowth (time between steps 8.9 and 8.11).

Note A: An alternative to step 8.11 involves the collection and counting of the strontium carbonate precipitate on a tared stainless steel planchet. For this, the approach is as follows:

1. Slurry the strontium carbonate precipitate with a few ml water and transfer quantitatively to a tared stainless steel planchet. Dry under infrared lamp.
2. Cool, weigh, and beta count (within 2 hours).

Note B: The calculation of the total strontium activity D , in the sample at this point in time can be made as follows:

$$D = \frac{C}{2.22 * EVR}$$

where:

C = net count rate, cpm,
 E = counter efficiency, for strontium-90
 V = liters of sample used,

R = fractional chemical yield, and
2.22 = conversion factor from dpm/pCi.

Strontium-90 (By Yttrium-90)

- 8.12 After counting the strontium carbonate for strontium-89 and strontium-90 activity, store the filter or the planchet for a measured period of ingrowth, then proceed with the following steps for yttrium-90 separation. A 2-week or longer ingrowth period is recommended for samples with very low strontium-90 activity. Step 8.9 was the beginning of this ingrowth period.
- 8.13 Undo the mylar covering from the nylon ring and disc, and transfer the filter to a small funnel which has been placed to drain into a 50 ml centrifuge tube. Dissolve the strontium precipitate by carefully wetting the filter with 5 ml of 6N HNO₃. Wet the filter with 2.0 ml yttrium carrier. Rinse the strontium and yttrium into the centrifuge tube by washing the filter with four 5-ml portions of 1N HNO₃. Remove the funnel from the centrifuge tube, discard the filter, and add 1 drop of wetting agent solution to the centrifuge tube. Swirl the tube to mix the contents thoroughly.
- Note: In the case of the stainless steel planchet:
1. After the period for yttrium-90 ingrowth, slurry the precipitate on the planchet with 2 ml water and transfer to a centrifuge tube with the aid of a rubber policeman. To make the transfer quantitative, wash the residue from the planchet with a small amount of 1N HNO₃. Dissolve the precipitate in the tube with sufficient 1N HNO₃, and dilute with water to 10 ml.
 2. Add 2.0 ml yttrium carrier and stir.
 3. Boil to expel dissolved carbon dioxide; cool to room temperature.
- 8.14 Precipitate the yttrium as hydroxide by adding 10 ml 15N NH₄OH to the centrifuge tube, stirring and heating for 10 minutes in a hot water bath. Cool, centrifuge and decant supernate into a 100-ml beaker. Note time of last precipitation; this is the end of yttrium-90 ingrowth and the beginning of yttrium-90 decay.
- 8.15 Dissolve precipitate in 1 ml 1N HNO₃ and dilute with water to 10 ml.
- 8.16 Reprecipitate yttrium by dropwise addition of 15N NH₄OH.
- 8.17 Centrifuge and combine supernate with solution in the 100-ml beaker (step 8.14).
- 8.18 Repeat steps 8.15, and 8.16. Save the combined supernates in the beaker for strontium gravimetric yield determination, step 8.22.

Note: Steps 8.22 to 8.25 are a repeat of the strontium carbonate precipitation to determine chemical yield after the yttrium has been removed.

- 8.19 Dissolve the precipitate in 15 ml water containing 2 ml 6N HCl. Precipitate the yttrium as oxalate by adding 20 ml saturated $(\text{NH}_4)_2\text{C}_2\text{O}_4$ and heating for 30 minutes in a hot water bath (near boiling). Cool in an ice bath and then filter the yttrium oxalate onto a Whatman #42 filter (4.25 cm diameter). Wash the precipitate with three 5-ml portions of water, then with three 5-ml portions of acetone. Air dry the filter for about 1 hour.

Note: A pH of 1.7-1.9 in the solution from which yttrium oxalate is being precipitated is necessary to get a uniform $9\text{H}_2\text{O}$ hydrate precipitate. This is necessary if the analyst prefers and is going to weigh the yttrium oxalate for chemical yield. Also, the analyst may then prefer to use a tared glass fiber filter instead of a Whatman #42 paper filter. The filter plus oxalate precipitate is weighed to determine chemical yield (recovery). See note following 9.1 for calculations. If this procedure is followed, step 8.21 can be eliminated.

- 8.20 Mount filter on a plastic ring and disc, and count for yttrium-90 activity. Record the time of the counting for decay correction (time between 8.14 and count time).
- 8.21 Undo the mylar covering, and transfer the filter to a tared crucible. Ignite at 800°C for 1 hour in a muffle furnace to convert the oxalate to the oxide. Cool and weigh the crucible. Determine the yttrium recovery (see Section 9.1).
- 8.22 Warm the combined supernates from step 8.14, add 5 ml 2N Na_2CO_3 , and digest for 10 minutes. Cool, centrifuge, and discard supernate.
- 8.23 Wash the SrCO_3 with 15 ml water and discard wash solution.
- 8.24 Slurry with a few ml water and transfer quantitatively to a tared stainless-steel planchet. Dry under infrared lamps.
- 8.25 Cool and weigh the planchet. Determine the strontium recovery (see Section 9.1).

9. Calculations

- 9.1 Chemical yields for strontium and yttrium

$$a = \text{Yield factor for Sr} = \frac{\text{mg SrCO}_3 \text{ recovered}}{\text{mg Sr carrier added (as carbonate)}}$$

20.0 mg of strontium is equivalent to 33.7 mg SrCO_3

$$b = \text{Yield factor for Y} = \frac{\text{mg Y}_2\text{O}_3 \text{ Recovered}}{\text{mg Y carrier added (as oxide)}}$$

20.0 mg of yttrium is equivalent to 25.4 mg of Y_2O_3 .

Note: If chemical yield is to be determined from the yttrium oxalate precipitate the following calculations are used.

$$b = \text{yield factor for Y} = \frac{\text{mg Y}_2(\text{C}_2\text{O}_4)_3 \text{CH}_2\text{O}}{\text{mg Y carrier added (as oxalate)}}$$

20.0 mg of yttrium is equivalent to 67.9 mg of yttrium oxalate, $\text{Y}_2(\text{C}_2\text{O}_4)_3 \cdot \text{H}_2\text{O}$.

9.2 Calculations for Activities at Equilibrium Conditions:

Indicated cpm values are net cpm (reagent blank, including background, subtracted).

$$\text{Y}^{90} \text{ dpm} = \frac{\text{Y}^{90} \text{ dpm}}{\text{abefi}} = \text{Sr}^{90} \text{ dpm}$$

$$\text{Sr}^{90} \text{ cpm} = \text{Y}^{90} \text{ dpm} * c$$

Total $\text{Sr}^{89,90}$ cpm = total cpm (SrCO_3) - ingrown Y^{90} cpm

$$\text{Ingrown Y}^{90} \text{ cpm} = \text{Sr}^{90} \text{ cpm} * e * g$$

$$\text{Sr}^{89} \text{ cpm} = \frac{\text{total Sr cpm}}{a} - \text{Sr}^{90} \text{ cpm} - \text{Sr}^{90} \text{ dpm} * e * g$$

$$\text{Sr}^{89} \text{ dpm} = \frac{\text{Sr}^{89} \text{ cpm}}{d}$$

$$\text{Sr}^{90} \text{ pCi/liter} = \frac{\text{Sr}^{90} \text{ dpm}}{2.22 * V}$$

$$\text{Sr}^{89} \text{ pCi/liter} = \frac{\text{Sr}^{89} \text{ dpm}}{2.22 * V}$$

where:

- c = strontium-90 counting efficiency,
- d = strontium-89 counting efficiency,
- e = yttrium-90 counting efficiency,

| | | |
|------|---|--|
| f | = | yttrium-90 decay factor, |
| g | = | yttrium-90 ingrowth factor, for unwanted yttrium-90 in total strontium-89, strontium-90 count, |
| h | = | strontium-89 decay factor, |
| V | = | volume of sample analyzed, in liters, |
| i | = | yttrium-90 ingrowth factor for strontium-90 by yttrium-90 determination, and |
| 2.22 | = | conversion factor from dpm/pCi. |

Error associated with the results of the analysis should also be reported. See Section 10 for error and statistical calculations, for yttrium-90 decay and ingrowth factors, and for strontium-89 decay factors.

10. Calculation Factors

10.1 Error and Statistical Calculations - Because of the random nature of radioactivity disintegrations, there is an error associated with any measured count of these disintegrations. The variability of any measurement is indicated by the standard deviation. The standard deviation in the counting rate, (R), is determined by the following equation:

$$\sigma (R) = \left[\frac{R_0}{t_1} + \frac{B}{t_2} \right]^{\frac{1}{2}}$$

where: R_0 = Gross count rate
 T_1 = counting time for the gross count
 B = background count rate
 T_2 = counting time for the background count
 Let $\sigma(R1) = D_1$ = standard deviation for the count of total strontium-89 and strontium-90 (from $SrCO_3$ precipitate, which includes the unwanted ingrown yttrium-90).

Let $\sigma(R2) = D_2$ = Standard deviation for the yttrium-90 count for the strontium-90 determination.

The counting errors, E, for a given sample for the strontium-89 and strontium-90 determinations expressed in pCi/liter are shown as follows:

$$\text{For } ^{90}\text{Sr}, E = \frac{1.96D_2 * 1000}{2.22 * abefiV}$$

$$\text{For } ^{89}\text{Sr}, E = \frac{1.96 * 1000}{2.22 * adV} \left[(D_1)^2 + \left(\frac{c + e * g}{befi} \right)^2 (D_2)^2 \right]^{\frac{1}{2}}$$

where:

- 1.96 = 95% confidence factor
- 1000 = ml/liter
- 2.22 = conversion factor from disintegrations/minute to picocuries
- a = strontium recovery factor
- b = yttrium recovery factor
- c = strontium-90 counting efficiency
- d = strontium-89 counting efficiency
- e = yttrium-90 counting efficiency
- f = yttrium-90 decay factor
- g = ingrowth factor for unwanted yttrium-90 in total radiostrontium count
- i = ingrowth factor for yttrium-90 for strontium-90 determination

These were derived by applying propagation of error theory to the expressions in Section 9.3.

The standard deviations of a number of experimental analyses or observations is determined by:

$$S = \left[\sum_{i=1}^n \frac{(x_i - \bar{x})^2}{n - 1} \right]^{\frac{1}{2}}$$

where:

- x_i = activity (pCi/liter) of a given sample
- \bar{x} = mean activity (pCi/liter) of a series of analyses
- n = the number of replicate analyses

10.2 Yttrium-90 Decay and Ingrowth Factors (0-71 Hours)

| t(hr) | $e^{-\lambda t}$ | $1-e^{-\lambda t}$ | t(hr) | $e^{-\lambda t}$ | $1-e^{-\lambda t}$ | t(hr) | $e^{-\lambda t}$ | $1-e^{-\lambda t}$ |
|-------|------------------|--------------------|-------|------------------|--------------------|-------|------------------|--------------------|
| 0.0 | 1.0000 | .0000 | 22.5 | .7843 | .2157 | 45.0 | .6151 | .3849 |
| 0.5 | .9940 | .0054 | 23.0 | .7801 | .2199 | 45.5 | .6118 | .3882 |
| 1.0 | .9893 | .0107 | 23.5 | .7759 | .2241 | 46.0 | .6085 | .3915 |
| 1.5 | .9839 | .0161 | 24.0 | .7717 | .2283 | 46.5 | .6053 | .3947 |
| 2.0 | .9786 | .0214 | 24.5 | .7676 | .2324 | 47.0 | .6020 | .3980 |
| 2.5 | .9734 | .0266 | 25.0 | .7634 | .2366 | 47.5 | .5988 | .4012 |
| 3.0 | .9681 | .0319 | 25.5 | .7593 | .2407 | 48.0 | .5955 | .4045 |
| 3.5 | .9629 | .0371 | 26.0 | .7552 | .2448 | 48.5 | .5923 | .4077 |
| 4.0 | .9577 | .0423 | 26.5 | .7512 | .2488 | 49.0 | .5891 | .4109 |
| 4.5 | .9526 | .0474 | 27.0 | .7471 | .2529 | 49.5 | .5860 | .4140 |
| 5.0 | .9474 | .0526 | 27.5 | .7431 | .2569 | 50.0 | .5828 | .4172 |
| 5.5 | .9423 | .0577 | 28.0 | .7391 | .2609 | 50.5 | .5797 | .4203 |
| 6.0 | .9373 | .0627 | 28.5 | .7351 | .2549 | 51.0 | .5766 | .4234 |
| 6.5 | .9322 | .0678 | 29.0 | .7311 | .2689 | 51.5 | .5735 | .4265 |
| 7.0 | .9272 | .0728 | 29.5 | .7272 | .2728 | 52.0 | .5704 | .4296 |
| 7.5 | .9222 | .0778 | 30.0 | .7233 | .2767 | 52.5 | .5673 | .4327 |
| 8.0 | .9172 | .0828 | 30.5 | .7194 | .2806 | 53.0 | .5642 | .4358 |
| 8.5 | .9123 | .0877 | 31.0 | .7155 | .2845 | 53.5 | .5612 | .4388 |
| 9.0 | .9074 | .0926 | 31.5 | .7117 | .2883 | 54.0 | .5582 | .4418 |
| 9.5 | .9025 | .0975 | 32.0 | .7078 | .2922 | 54.5 | .5552 | .4448 |
| 10.0 | .8976 | .1024 | 32.5 | .7040 | .2960 | 55.0 | .5522 | .4478 |
| 10.5 | .8928 | .1072 | 33.0 | .7002 | .2998 | 55.5 | .5492 | .4508 |
| 11.0 | .8880 | .1120 | 33.5 | .6965 | .3035 | 56.0 | .5462 | .4538 |
| 11.5 | .8832 | .1168 | 34.0 | .6927 | .3073 | 56.5 | .5433 | .4567 |
| 12.0 | .8785 | .1215 | 34.5 | .6890 | .3110 | 57.0 | .5404 | .4596 |
| 12.5 | .8737 | .1263 | 35.0 | .6853 | .3147 | 57.5 | .5375 | .4625 |
| 13.0 | .8690 | .1310 | 35.5 | .6816 | .3184 | 58.0 | .5346 | .4654 |
| 13.5 | .8644 | .1356 | 36.0 | .6779 | .3221 | 58.5 | .5317 | .4683 |
| 14.0 | .8597 | .1403 | 36.5 | .6743 | .3257 | 59.0 | .5288 | .4712 |
| 14.5 | .8551 | .1449 | 37.0 | .6706 | .3294 | 59.5 | .5260 | .4740 |
| 15.0 | .8505 | .1495 | 37.5 | .6670 | .3330 | 60.0 | .5232 | .4768 |
| 15.5 | .8459 | .1541 | 38.0 | .6634 | .3366 | 60.5 | .5203 | .4797 |
| 16.0 | .8413 | .1587 | 38.5 | .6599 | .3401 | 61.0 | .5175 | .4825 |
| 16.5 | .8468 | .1632 | 39.0 | .6563 | .3437 | 61.5 | .5148 | .4852 |
| 17.0 | .8323 | .1677 | 39.5 | .6528 | .3472 | 62.0 | .5092 | .4880 |
| 17.5 | .8278 | .1722 | 40.0 | .6493 | .3507 | 62.5 | .5092 | .4908 |
| 18.0 | .8234 | .1766 | 40.5 | .6458 | .3542 | 63.0 | .5065 | .4935 |
| 18.5 | .8189 | .1811 | 41.0 | .6423 | .3577 | 64.0 | .5010 | .4990 |
| 19.0 | .8145 | .1855 | 41.5 | .6388 | .3612 | 65.0 | .4957 | .5043 |
| 19.5 | .8101 | .1899 | 42.0 | .6354 | .3646 | 66.0 | .4903 | .5097 |
| 20.0 | .8058 | .1942 | 42.5 | .6320 | .3680 | 67.0 | .4851 | .5149 |
| 20.5 | .8014 | .1986 | 43.0 | .6286 | .3714 | 68.0 | .4799 | .5201 |
| 21.0 | .7971 | .2029 | 43.5 | .6252 | .3748 | 69.0 | .4747 | .5253 |
| 21.5 | .7928 | .2072 | 44.0 | .6219 | .3781 | 70.0 | .4696 | .5304 |
| 22.0 | .7885 | .2115 | 44.5 | .6185 | .3815 | 71.0 | .4646 | .5354 |

10.3 Yttrium-90 Ingrowth Factors (0-27 days)

| t(days) | 1-e- λ t | t(days) | 1-e- λ t | t(days) | 1-e- λ t |
|---------|------------------|---------|------------------|---------|------------------|
| 0.00 | .0000 | 9.00 | .9029 | 18.00 | .9906 |
| 0.25 | .0627 | 9.25 | .9090 | 18.25 | .9912 |
| 0.50 | .1215 | 9.50 | .9147 | 18.50 | .9917 |
| 0.75 | .1766 | 9.75 | .9201 | 18.75 | .9922 |
| 1.00 | .2283 | 10.00 | .9251 | 19.00 | .9927 |
| 1.25 | .2767 | 10.25 | .9298 | 19.25 | .9932 |
| 1.50 | .3221 | 10.50 | .9342 | 19.50 | .9936 |
| 1.75 | .3646 | 10.75 | .9384 | 19.75 | .9940 |
| 2.00 | .4045 | 11.00 | .9422 | 20.00 | .9944 |
| 2.25 | .4418 | 11.25 | .9458 | 20.25 | .9948 |
| 2.50 | .4768 | 11.50 | .9492 | 20.50 | .9951 |
| 2.75 | .5097 | 11.75 | .9524 | 20.75 | .9954 |
| 3.00 | .5404 | 12.00 | .9554 | 21.00 | .9957 |
| 3.25 | .5692 | 12.25 | .9582 | 21.25 | .9959 |
| 3.50 | .5963 | 12.50 | .9608 | 21.50 | .9962 |
| 3.75 | .6216 | 12.75 | .9633 | 21.75 | .9964 |
| 4.00 | .6453 | 13.00 | .9656 | 22.00 | .9967 |
| 4.25 | .6676 | 13.25 | .9678 | 22.25 | .9969 |
| 4.50 | .6884 | 13.50 | .9697 | 22.50 | .9971 |
| 4.75 | .7080 | 13.75 | .9716 | 22.75 | .9973 |
| 5.00 | .7263 | 14.00 | .9734 | 23.00 | .9974 |
| 5.25 | .7435 | 14.25 | .9751 | 23.25 | .9976 |
| 5.50 | .7596 | 14.50 | .9766 | 23.50 | .9977 |
| 5.75 | .7746 | 14.75 | .9781 | 23.75 | .9979 |
| 6.00 | .7888 | 15.00 | .9795 | 24.00 | .9980 |
| 6.25 | .8020 | 15.25 | .9808 | 24.25 | .9981 |
| 6.50 | .8145 | 15.50 | .9820 | 24.50 | .9982 |
| 6.75 | .8261 | 15.75 | .9831 | 24.75 | .9984 |
| 7.00 | .8370 | 16.00 | .9842 | 25.00 | .9985 |
| 7.25 | .8472 | 16.25 | .9852 | 25.25 | .9986 |
| 7.50 | .8568 | 16.50 | .9861 | 25.50 | .9987 |
| 7.75 | .8658 | 16.75 | .9870 | 25.75 | .9987 |
| 8.00 | .8742 | 17.00 | .9878 | 26.00 | .9988 |
| 8.25 | .8820 | 17.25 | .9886 | 26.25 | .9989 |
| 8.50 | .8896 | 17.50 | .9893 | 26.50 | .9990 |
| 8.75 | .8864 | 17.75 | .9900 | 26.75 | .9990 |
| | | | | 27.00 | .9991 |

10.4 Strontium-89 Decay Factors (0-59.5 days) ($t_{1/2} = 51$ days)

| t(days) | $e^{-\lambda t}$ | t(days) | $e^{-\lambda t}$ | t(days) | $e^{-\lambda t}$ |
|---------|------------------|---------|------------------|---------|------------------|
| 0.0 | 1.0000 | 20.0 | .7620 | 40.0 | .5808 |
| 0.5 | .9932 | 20.5 | .7569 | 40.5 | .5769 |
| 1.0 | .9865 | 21.0 | .7518 | 41.0 | .5730 |
| 1.5 | .9798 | 21.5 | .7568 | 41.5 | .5690 |
| 2.0 | .9732 | 22.0 | .7416 | 42.0 | .5652 |
| 2.5 | .9668 | 22.5 | .7366 | 42.5 | .5613 |
| 3.0 | .9601 | 23.0 | .7317 | 43.0 | .5575 |
| 3.5 | .9536 | 23.5 | .7267 | 43.5 | .5539 |
| 4.0 | .9471 | 24.0 | .7218 | 44.0 | .5500 |
| 4.5 | .9407 | 24.5 | .7169 | 44.5 | .5462 |
| 5.0 | .9344 | 25.0 | .7120 | 45.0 | .5427 |
| 5.5 | .9280 | 25.5 | .7072 | 45.5 | .5380 |
| 6.0 | .9217 | 26.0 | .7023 | 46.0 | .5352 |
| 6.5 | .9155 | 26.5 | .6977 | 46.5 | .5318 |
| 7.0 | .9093 | 27.0 | .6930 | 47.0 | .5280 |
| 7.5 | .9031 | 27.5 | .6882 | 47.5 | .5245 |
| 8.0 | .8970 | 28.0 | .6836 | 48.0 | .5210 |
| 8.5 | .8909 | 28.5 | .6790 | 48.5 | .5175 |
| 9.0 | .8849 | 29.0 | .6742 | 49.0 | .5140 |
| 9.5 | .8789 | 29.5 | .6699 | 49.5 | .5105 |
| 10.0 | .8729 | 30.0 | .6651 | 50.0 | .5070 |
| 10.5 | .8670 | 30.5 | .6608 | 50.5 | .5035 |
| 11.0 | .8612 | 31.0 | .6562 | 51.0 | .5000 |
| 11.5 | .8553 | 31.5 | .6519 | 51.5 | .4967 |
| 12.0 | .8495 | 32.0 | .6473 | 52.0 | .4933 |
| 12.5 | .8438 | 32.5 | .6430 | 52.5 | .4900 |
| 13.0 | .8381 | 33.0 | .6388 | 53.0 | .4868 |
| 13.5 | .8324 | 33.5 | .6342 | 53.5 | .4834 |
| 14.0 | .8268 | 34.0 | .6300 | 54.0 | .4801 |
| 14.5 | .8212 | 34.5 | .6259 | 54.5 | .4769 |
| 15.0 | .8156 | 35.0 | .6215 | 55.0 | .4734 |
| 15.5 | .8101 | 35.5 | .6172 | 55.5 | .4702 |
| 16.0 | .8046 | 36.0 | .6131 | 56.0 | .4671 |
| 16.5 | .7991 | 36.5 | .6090 | 56.5 | .4640 |
| 17.0 | .7938 | 37.0 | .6050 | 57.0 | .4608 |
| 17.5 | .7883 | 37.5 | .6009 | 57.5 | .4578 |
| 18.0 | .7881 | 38.0 | .5968 | 58.0 | .4547 |
| 18.5 | .7778 | 38.5 | .5928 | 58.5 | .4513 |
| 19.0 | .7725 | 39.0 | .5888 | 59.0 | .4484 |
| 19.5 | .7672 | 39.5 | .5848 | 59.5 | .4454 |

11. Precision and Accuracy

- 11.1 In a single operator test of the method, two sets of five water samples containing known amounts of strontium-89, strontium-90, were analyzed for those radionuclides. The average recovery of added strontium-90 was 95 and 94 percent for the two sets of samples at a precision of 3 and 5 percent at the 95 percent confidence level.
- 11.2 In a collaborative test of the method with 13 laboratories participating, three samples containing known amounts of strontium-89 and strontium-90 were analyzed (Samples A, B and C).
- 11.3 The data of two laboratories for all three samples for strontium-90 were rejected because their scores in the ranked results of the laboratory averages were outside the acceptable range for 13 laboratories and 3 samples.
- 11.4 The coefficients of variation for the three samples ranged from 11.3% for 1000 pCi/l concentrations to 57% for 10 pCi/l concentrations.
- 11.5 The coefficients of variation for the combined within-laboratory precision for strontium-90 in the three samples ranged from 13.6% for 1000 pCi/l concentrations to 23% for 10 pCi/l concentrations.
- 11.6 The coefficients of variation for the precision of the method between laboratories for strontium-89 in the three samples ranged from 20% for 1000 pCi/l concentrations to 43% for 10 pCi/l concentrations.
- 11.7 The coefficients of variation for the precision of the method between laboratories for strontium-90 in the samples ranged from 25% for 1000 pCi/l concentrations to 44% for 10 pCi/l concentrations.
- 11.8 The coefficients of variation for the total error between laboratories based on a single analysis for strontium-89 in the three samples ranged from 23% for 1000 pCi/l concentrations to 71% for 10 pCi/l concentrations.
- 11.9 The coefficients of variation for the total error between laboratories based on a single analysis for strontium-90 in the three samples ranged from 17% for 1000 pCi/l concentrations to 46% for 10 pCi/l concentrations.
- 11.10 In the statistical test to detect method bias, no significant bias was shown in the analysis of the three samples for strontium-89. In the analysis of the three samples for strontium-90, sample C (1000 pCi/l concentration) showed a low bias but not seriously.
- 11.11 The strontium-89 analysis of samples A, B, and C deviated from the known values by the factor 1.49, 1.01, and 1.03, respectively. The strontium-90 analysis of

samples A, B, and C deviated from the known values by the factors 1.00, 0.899, and 0.820, respectively.

Bibliography

1. EPA Drinking Water Regulations, Federal Register, Vol. 41, No. Friday July 9, 1976, Part 141.
2. American Public Health Association. Standard Methods for the Examination of Water and Wastewater, 13th Edition, Washington, D.C. 1971.
3. American Society for Testing and Materials. 1976 Annual Book of ASTM Standards, Part 31, Philadelphia, PA.
4. Dixon, W.J., and F. J. Massey, Jr. Introduction to Statistical Analysis, McGraw-Hill, 1969.
5. HASL-300, Health and Safety Laboratory, ERD, 376 Hudson Street, New York, New York 10014.
6. Steiner, E.H., 1975. Planning and Analysis for Results of Collaborative Tests, Statistical Manual of the AOAC Association of Official Analytical Chemists, Washington, D. C.
7. Youden, W. J. 1969. Precision Measurements and Calibration. H. H. Ku, editor. National Bureau of Standards, Special Publication 300, Vol. 1.
8. Youden, W. J. 1975. Statistical Techniques for Collaborative Tests. Statistical Manual of the AOAC. Association of Official Analytical Chemists, Washington, D. C.
9. Douglas, G. S., ed. Radioassay Procedures for Environmental Samples, Environmental Health Series, USDHM Rept. 999-RH-27, National Center for Radiological Health, Rockville, Md. 10852 (Jan. 1967).
10. Hahn, R.B., and C. P. Straub. "Determination of Radioactive Strontium and Barium in Water." J. Am. Water Works Assoc. 47 (4) 335-340 (April 1955)
11. Goldin, A. S., R. J. Velten and G. W. Frishkorn, "Determination of Radioactive Strontium," Anal. Chem. 31, 1490 (1959).